**VOCO**

## LC Dam

**(GB)**

**light-curing gingival protection Instructions for use**

**(CE)**

**LC Dam** is a light-curing composite for covering the gingiva and neighbouring teeth during whitening therapy. It provides a barrier against highly concentrated whitening agents and protects the surrounding soft tissue and teeth that are not to be whitened.

**Fields of application:**

**LC Dam** covers and protects the gingiva and neighbouring teeth during whitening therapy carried out in the dental surgery (in-office whitening).

**Application:**

Before starting treatment, allow the material to reach ambient temperature and insert the cheek retractor. Create a relatively dry treatment area in the usual way, e.g. with cotton rolls and air syringe. Remove the protective cap and attach a cannula to the syringe.

Apply **LC Dam** onto the gingiva adjoining the tooth in sections of two to three teeth, light cure each section briefly with the polymerisation lamp and then completely light-cure each section. Conventional polymerisation devices are suitable for light-curing the material. The curing time is 20 seconds per light guide diameter for a lamp with a light output of 1000 mW/cm². Shorten the duration of exposure accordingly for lamps with a higher power density and extend the duration of exposure accordingly for lamps with a lower power density. A 4 - 6 mm wide strip of 1.5 to 2 mm thickness will suffice to ensure that the barrier provides effective protection. About 0.5 mm of the cervical dental hard tissue must also be covered to ensure effective protection of the gingiva. Remove any large excess amounts of **LC Dam** from the tooth enamel with a suitable hand instrument (e.g. probe, scaler or curette) to prevent impairment of the whitening effect of **Perfect Bleach Office**.

**LC Dam** must be applied up to and including the first tooth not being whitened. The interdental spaces must also be sealed with **LC Dam** to protect the interdental papillae from oxidation. Reseal the syringe with the protective cap after application.

Whitening treatment can begin immediately after **LC Dam** has been carefully light-cured (see above).

**LC Dam** can be easily removed with an instrument (see above) on completion of the whitening therapy and the removal of all residual whitening gel. The material can normally be completely removed in a few large pieces. Remove any remaining material from the interdental areas with a probe or dental floss. Allow patients to rinse their mouths thoroughly with water afterwards.

**Information, indications, precautions, contraindications:** **LC Dam** is contraindicated for patients with known hypersensitivities to product ingredients (methacrylate and BHT). If irritation of an allergic reaction should occur, stop the treatment immediately, remove **LC Dam** and rinse the treatment area with plenty of water. Refer patients to a doctor, if necessary. Phenolic substances, especially products containing eugenol and thymol, lead to curing failures and are to be avoided.

**Storage:**

Store **LC Dam** at a temperature of between 2°C - 23°C (refrigerator or ambient temperature). Do not use **LC Dam** after its expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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**VOCO GmbH**